

Included in package:



1 Fabrication Dummy



CD103FD

Fabrication Dummy

For CD103 Air-Lock or CD117 Easy Off Lock
CD103D Deep Air-Lock or CD117D Easy Off Deep lock



Manufactured by:



419 N. Curtis Rd., ID 83706 USA
(208) 429-0026
www.coyote.us

2-Year warranty against manufacturer defects, excessive wear or breakage

External Prosthetic Components



EC REP

Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta



CD103FD.revA.04162021

Included in package:



1 Fabrication Dummy



CD103FD

Fabrication Dummy

For CD103 Air-Lock or CD117 Easy Off Lock
CD103D Deep Air-Lock or CD117D Easy Off Deep lock



Manufactured by:



419 N. Curtis Rd., ID 83706 USA
(208) 429-0026
www.coyote.us

2-Year warranty against manufacturer defects, excessive wear or breakage

External Prosthetic Components



EC REP

Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta



CD103FD.revA.04162021

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103FD



1. Put CD103FD fabrication dummy over anchor.



2. Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



3. Pressing in around fabrication dummy will help remove any air pockets.



4. If needed use parachute cord or rope to pull plastic tight to dummies edge.



5. Let plastic cool before cutting fabrication dummy free.



6. Cut tooling free with cast saw, knife or utility knife.



7. Remove tooling from anchor.



8. Sand flexible inner down so anchor is almost 100% exposed



9. Pull a PVA bag over the anchor and cast. Glue your lock to anchor and start laminating.

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103FD



1. Put CD103FD fabrication dummy over anchor.



2. Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



3. Pressing in around fabrication dummy will help remove any air pockets.



4. If needed use parachute cord or rope to pull plastic tight to dummies edge.



5. Let plastic cool before cutting fabrication dummy free.



6. Cut tooling free with cast saw, knife or utility knife.



7. Remove tooling from anchor.



8. Sand flexible inner down so anchor is almost 100% exposed



9. Pull a PVA bag over the anchor and cast. Glue your lock to anchor and start laminating.